

FAST, ACCURATE AND
RELIABLE **LACTATE**
MEASUREMENT



Easy to use

Fast and accurate
results

Performance
measurement

Practical and
reliable

EKF
DIAGNOSTICS

SCIENCE MADE SIMPLE



<p>EASY TO USE</p>	<ul style="list-style-type: none"> • Insert/remove sensor to turn on/off • Jog dial navigation • Just 0.5µl of capillary blood required • Pre-calibrated sensors • Device calibrated by simple manual coding • Lactate measurement, date/time, mode, temperature and memory ID on a single display • Automatic self test
<p>FAST AND ACCURATE RESULTS</p>	<ul style="list-style-type: none"> • Enzymatic amperometric detection method • Results within 10 seconds • Measuring range: 0.5 – 25mmol/L • Imprecision: ≤ 3.4mmol/L; ± 0.2mmol/L (> 3.4mmol/L : ±3%) • Test solutions available for function control
<p>PERFORMANCE MEASUREMENT</p>	<ul style="list-style-type: none"> • Single and step-test measurements (resting/exercise/recreation) • Stopwatch and countdown timer functions • Lactate Scout Assistant performance management software available
<p>PRACTICAL AND RELIABLE</p>	<ul style="list-style-type: none"> • Pocket-size: 91mm (h) x 55mm (w) x 24mm (d) • Lightweight: 80g • 1,000 tests using just 2 x AAA batteries • Operating range: 5 – 45°C and 85% humidity • Stores up to 250 results • PC and Bluetooth connectivity

Simple, accurate and reliable lactate measurement in 3 easy steps



1

Remove sensor and place into analyser.



2

Prick finger and collect blood by touching sensor.



3

Result appears in 10 seconds.

Distributed by



Manufacturer

SensLab GmbH
 Bautzner Strasse 67
 D-04347 Leipzig
Tel +49 (0)341 2341841

Sales

EKF Diagnostics Holdings plc
 Avon House, 19 Stanwell Road,
 Penarth, Cardiff, UK, CF64 2EZ
Tel +44 (0)29 20 710570
www.ekfdiagnostics.com

Approved for the European Market



SCIENCE MADE SIMPLE

Disclaimer: EKF Diagnostics Holdings plc reserves the right to change the product specification without prior notice.

Revision 0.0-06/2011