

# McGRATH<sup>®</sup> MAC

## Video Laryngoscope

*Improve Your View*



Distributed by



**COVIDIEN**


# Improve your view with the McGRATH® MAC video laryngoscope

## THE ANTERIOR CAMERA IS DESIGNED TO

- Improve grade of view
- Provide better visual confirmation of tube placement
- Reduce tunnel vision experienced with direct laryngoscopy

## MINIMAL CHANGE IN INTUBATION TECHNIQUE

- Proprietary Macintosh type blade
- Used like a traditional laryngoscope





### INTEGRATED VERTICAL SCREEN

- Designed to reduce blind spots
- Keeps the focus on your patient

### TURN ON AND GO

- Single-button operation
- Immediately ready for use when you are

### DURABLE

- Drop tested up to two meters

### LONG-LASTING BATTERY

- 250 minutes
- Remaining time clearly displayed on screen

### EASY TO CLEAN

- (HLDi) High Level Disinfection-immersible
- Compatible with the latest vaporized hydrogen peroxide sterilization systems

### DISPOSABLE BLADES

- Packaged sterile
- Anti-fog coating

### DENTAL FRIENDLY

- Slim-blade design

**Designed for safety, simplicity and strength**

## Ordering information

### MCGRATH® MAC VIDEO LARYNGOSCOPE

SKU NO.	DESCRIPTION	SIZE	QUANTITY
300-000-000	Handle	NA	1
350-005-000	McGRATH® MAC Blade Size 3	3	50/box
350-013-000	McGRATH® MAC Blade Size 4	4	50/box
340-000-000	250-Minute Battery Pack	NA	1

### TECHNICAL SPECIFICATIONS

#### McGRATH® MAC Handle

Size:	180mm x 68mm x 110mm
Weight:	200g
Power:	Proprietary 3.6V Lithium Battery Pack (c. 250 minutes )
Light source:	High intensity LED
Display:	2.5" LCD color display
Camera:	CMOS
Materials:	Durable medical grade thermoplastics with reinforced structural alloy core. The device and packaging are latex free.

#### McGRATH® MAC Disposable Blades

Material:	Fog-free medical grade optical polymer
Packaging:	Packaged sterile for single use

For demonstrations and quotes, please contact your local representative. Details can be found at [solutions.covidien.com](http://solutions.covidien.com).



# McGRATH® MAC



Designed & manufactured  
by Aircraft Medical

"The McGRATH® MAC video laryngoscope complies with EN 60601-1 and EN 60601-1-2 safety standards. The CE mark indicates that it meets the requirements of European Council Directive 93/42/EEC concerning medical devices. The device is regulated in the USA under FDA Regulation Number 868.5540 and device listed under the name "McGRATH MAC".

"McGRATH" and "Aircraft" are registered trademarks of Aircraft Medical Limited. "CameraStick" is a trademark of Aircraft Medical Limited.

COVIDIEN, COVIDIEN with logo, Covidien logo and *positive results for life* are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company. ©2011 Covidien. All rights reserved.

11-AW-0248 AW23011

Distributed by  
 **COVIDIEN**

6135 GUNBARREL AVENUE  
BOULDER, CO  
80301  
800-635-5267

[WWW.COVIDIEN.COM](http://WWW.COVIDIEN.COM)